


Fit it snugly but comfortably over the groin

## WEAR & CARE INSTRUCTIONS




Start with the instructions for JOBSKIN®

working it up to the middle of your thigh

about 5cm at a time.

Knee Length garment (Fig. 1-7). Continue



Stand up. Gather the fabric in both hands

and pull it up to your waist.

Unit 1, Level 1, 100 Station Street Nunawading Victoria 3131, Australia AUT&F: 1800 772 577


NZT&F: 0800 765 627 E <u>customer.service@jobskin.com</u> W www.jobskin.com

## WEAR & CARE INSTRUCTIONS



13

Adjust your JOBSKIN® garment so that it is snug at your midsection. This will keep it from sliding down. Pull the midsection to its proper height.



14

Slip your hands, palm out, into the back and lift the fabric up and away from your skin. Let it fall back to fit snugly over buttocks.



15

Check whether the garment is left or right.
Open the zipper completely.



19

(This garment may be difficult to put on alone. You may need help). Follow the directions for JOBSKIN® zippered garment (Fig. 15-18)



20

Stand up and continue to close zipper holding the zippered area together 10-15cm above the zipper tab. Adjust the leg along the inside over your groin and buttocks



2

From this point, use the instructions for Waist High garment without zipper (Fig. 10, 11, 12) to ease your JOBSKIN® garment comfortably into place. Adjust garment if necessary.



25

Optional: Insert pad at bend of elbow for added comfort. Fold sleeve down to elbow. Tuck pad in at elbow. Smooth the sleeve up over the pad and continue to ease it up the arm.



26

Straps should be securely fastened.



27

(Assistance may be necessary). Open zipper fully. Place your fingers inside the sleeve at the "closed" end and stretch it as wide as possible. Slip the sleeve over the arm so that the black dot lies at the elbow. Refer to Fig. 25 for correct placement of optional pad.



3

Adjust the fingers by carefully pulling on the seams until the glove fits snugly into the web spaces between the fingers. Be sure the fabric is smoothed out evenly over the fingers and hand. Hold the hand up to the light. You should not be able to see light through the web space.



32

Close the zipper, if there is one.



33

Open Velcro® closure. Hold chin in place with one hand, if you have assistance, or press chin against upper sternum. Gather fabric and stretch over head.




16

Point your toes and slip your foot into the foot area. Place your thumb between the zipper and your leg to protect your skin.



17

Hold the zipper together with one hand approximately 10-15cm above the zipper tab.



18

Zip it closed, keeping the elastic backing between your skin and the zipper. Adjust any wrinkles or uncomfortable spots by gently pulling the fabric until the pressure is even.



22

Follow the steps for putting on garments with zippers. (Fig. 15) Ease your foot through the entire length of the garment. Refer to Fig. 10, 11, 12 to ease your JOBSKIN® garment comfortably into place.



Slip your hand through the entire sleeve keeping the seam in line with your little finger. Ease the sleeve up your arm. The black dot indicates where your elbow should be. It may be helpful to brace your hand against a wall.



2

If a gauntlet (fingerless glove) is attached, seam should run from the palm of your hand then curve over your elbow to your upper arm.



28

Close the zipper, holding the zipper together with one hand 10-15cm below the zipper tab. Be sure to keep the elastic backing between the skin and zipper. Continue down the arm to the wrist.




29

Method 1: Turn glove inside out. Fold tips of fingers inside. Place tips of glove on fingers. Work fabric down until thumb tip is covered. Turn hand of glove right side out over fingers. Method 1 is easier if your glove has a zipper.



30

Method 2: Slip your hand through wrist or glove and into fingers.



34

Adjust the facial features. The face mask must have the fabric distributed to fit properly. Adjustments may be made by slipping fingers through openings and placing fabric where desired.



35

Method 1: Follow directions for masks, (Fig. 33 & 34). Method 2: Close the neck portion around patient's neck. Pull head strap and back elastic over head. Refer to (Fig. 34) to make adjustments.



36

Always apply the garment on the more severe arm first (Refer to Fig. 23), when sleeves are fitted. Close the Velcro® closure or zipper. Adjust neck edge and waist for comfort.