




WASHING & CARE INSTRUCTIONS

JOBSKIN® garments are specifically prescribed and designed for each person. In order to keep the garments in good condition and to gain as much benefit as possible from the garments, please observe the following suggestions for their care:




Wash the garment **daily**. Hand or machine wash (using a laundry bag) in warm water. Separate light from dark colours. Wash bright colours separately before use. **Keep Velcro® closures closed when washing or storing garment.**



Do not tumble dry. Do not dry in the sun. Do not use a heater or dryer as this may cause deterioration of the fabric.



Do not use chlorine bleach. If the garment is worn in a chlorinated swimming pool or in salt water, rinse well before drying.



Do not iron.



Do not dryclean.



Water-based or silicone-based moisturisers may be used. Massage well into the skin before donning garment. Avoid wax, paraffin or petroleum-based moisturisers.



Use a mild liquid detergent and rinse well in warm water to remove all traces of detergent. Do not use high phosphate detergents.



Gently squeeze excess moisture out of the garment and roll in a towel to further remove moisture. Dry flat or drip dry in the shade.

Closely follow the application instructions which have been provided. Your therapist will demonstrate the most suitable way to don the garments. It is important to maintain regular contact with your therapist to ensure that the garment is fitting well. Small wounds should be covered with a dressing for protection, and to prevent soiling the garment.