

# READY-TO-WEAR GLOVES



Note : Glove fingers can be trimmed to length - no finishing required.

# Wear and Care Instructions

AUSTRALIAN MADE


Your therapist will demonstrate the most suitable way to don the garments. Closely follow the application instructions which have been provided. It is important to maintain regular contact with your therapist to ensure that the garment is fitting well. Two options for applying gloves are shown below.

## Sizing Information (Gloves) Available in Biege and Black.

| Size        | Wrist           | Palm            |
|-------------|-----------------|-----------------|
| Small       | 12.8cm - 15.3cm | 16.2cm - 19.3cm |
| Medium      | 15.3cm - 17.8cm | 19.3cm - 22.3cm |
| Large       | 17.8cm - 20.3cm | 22.3cm - 24.8cm |
| Extra Large | 20.3cm - 22.8cm | 24.8cm - 28.3cm |
| 2XL         | ????cm - ???cm  | ???cm - ???cm   |



Circumferential measurements are taken at the widest section of the hand (above the thumb eminence) and at the wrist crease (distal styloid process)



**Method 1 (preferred):** Turn glove inside out. Fold tips of fingers inside. Place tips of glove on fingers. Work fabric down until thumb tip is covered. Turn hand of glove right side out over fingers.



**Method 2:** Slip your hand through wrist or glove and into fingers.



Adjust the fingers by carefully pulling on the seams until the glove fits snugly into the web spaces between the fingers. Be sure the fabric is smoothed out evenly over the fingers and hand. Hold the hand up to the light. You should not be able to see light through the web space.



Wash the garment **daily**. Hand or machine wash (inside a laundry bag) in warm water. Separate light from dark colours. Keep Velcro® closures closed when washing or storing garment.



Do not tumble dry. Do not dry in the sun. Do not use a heater or dryer as this may cause deterioration of the fabric.




Do not use chlorine bleach. If the garment is worn in a chlorinated swimming pool or in salt water, rinse well before drying.



Do not iron.



Do not dryclean.



Water-based or silicone-based moisturisers may be used. Massage well into the skin before donning garment. Avoid wax, paraffin or petroleum-based moisturisers.



Use a mild liquid detergent and rinse well in warm water to remove all traces of detergent. Do not use high phosphate detergents.



Gently squeeze excess moisture out of the garment and roll in a towel to further remove moisture. Dry flat or drip dry in the shade.