

Patient details (Please print in black with CAPITAL letters)

Working together for better solutions

+61 3 9915 8000
customer.service@jobskin.com
Jobskin.com

Last name	First name	Date
-----------	------------	------



Measurements											
Pleat				Pleat				Toe end			
Left		Right		Left		Right		Left		Right	
mm		mm		mm		mm		mm		mm	
								-6			
								-4½			
								-3			
								-1½			
								Heel 0			
								1½			
								3			
								4½			
								6			
								7½			
								9			
								10½			
								12			
								13½			
								15			
Pleat				Knee End				Pleat			

Note: Seamless sole standard for 10 years and under
Note: Only pleat tapes at ends ie: proximal and/or distal end of garment

Personalisation		
Fabric Selection <input type="checkbox"/> Conformtex <input type="checkbox"/> Powernet <input type="checkbox"/> Powertex	Colour Selection Refer to Jobskin swatch or website for fabric and thread colour choices	
Fabric Colour	Thread Colour	Qty
1	1	
2	2	
3	3	

Comments
Photos assist in garment design

Garment styles			
Code	L	R	Description
0015			Calf Sleeve
0101			Knee Length
0105			Anklet (Garment length up to + 4 1/2 tape)
0041			Interdigital Web Spacer (to be worn over Foot Glove 0538)
1198B			Toe Piece <input type="checkbox"/> Open Tips <input type="checkbox"/> Closed Tips
0539			Thong Toe Anklet <input type="checkbox"/> Open Tips <input type="checkbox"/> Closed Tips
0539			Thong Toe Knee Length <input type="checkbox"/> Open Tips <input type="checkbox"/> Closed Tips
0538			Foot Glove <input type="checkbox"/> Open Tips <input type="checkbox"/> Closed Tips
0538			Individual Toe, Knee Length <input type="checkbox"/> Open Tips <input type="checkbox"/> Closed Tips

Options and modifications			
Code	L	R	Description
1169			Slant Inserts
0000			5cm Elastic
SBE			Silicon Backed Elastic <input type="checkbox"/> 2.5cm <input type="checkbox"/> 5cm
1190			Anterior Ankle Insert (powernet & lining)
1201			Seamless Sole (Standard up to the age of 10)
1183			Suede Reinforcement <input type="checkbox"/> Heel <input type="checkbox"/> Ball <input type="checkbox"/> Full
1159			Self Enclosed Toe
1160			Soft Enclosed Toe
1191			Slingback (for toe piece)
1164			Zippers (anklet)
1164			Zipper (knee-high)
ZCT			Zipper Velcro® Closure Tab (top of zipper)
1163B			Zipper Fastening Aids (along length of zip)
1172A			Zipper Loop <input type="checkbox"/> Leather <input type="checkbox"/> Ribbon
0000			Zipper Lining <input type="checkbox"/> Satin <input type="checkbox"/> Elastic
0027			Padded Zipper Lining
1187			Lining-Heel
1189			Lining-Anterior Ankle
1188			Lining-Elsewhere (specify location) <input type="checkbox"/> Satin <input type="checkbox"/> Cotton
1195			Fully Lined Garment <input type="checkbox"/> Soft Satin <input type="checkbox"/> Cotton
			SilonTEX® lining - Indicate location and size on either foot outline, TBSA map or printed photo
0027			Padding - Soft (Indicate on TBSA Map) <input type="checkbox"/> Pocket <input type="checkbox"/> Attached
0027			Padding - Velfoam (Indicate on TBSA Map) <input type="checkbox"/> Pocket <input type="checkbox"/> Attached
1202			Motif

Leg zipper locations				
	Location		Full Length	Or Specify Length (cm)
	L	R		
Lateral				
Medial				
Other				

Knee High, Anklet, Toe Piece

Measurement Form - Page 2 of 3

Patient details (Please print in black with CAPITAL letters)

Last name

First name

Date



LEFT

- Open Tips
- Closed Tips


RIGHT

- Open Tips
- Closed Tips

White boxes are lengths and grey boxes are circumferences

TOE CIRCUMFERENCES

Toe circumference 'B' should be 'base of toe' & 'A' should be 'at point of finished length' (mm), or widest point of toe.



Note: FOOT OUTLINE (SEE: PG 3) MUST ACCOMPANY THIS FORM

Knee High, Anklet, Toe Piece

Please ensure this page prints to correct size BEFORE taking foot trace.




Measurement Form - Page 3 of 3

Patient details (Please print in black with CAPITAL letters)

Last name

First name

Date



Velcro is a registered trademark of VELCRO BVBA

Design options




Foot outline measuring tips

Centre the patient's foot on the foot outline page as marked. Ensure that the patient is standing or sitting with the foot flat on the ground. Mark the following clearly:


- Back of the heel
- Webspaces
- Toe lengths (if measuring a toepiece or glove)
- Where the garment will finish

Use the grey JOBSKIN leg tape for these garments.
Use the pink JOBSKIN tape for paediatrics.




Pleat instructions

When pleating the last tape, measure the distance between the pleated tape and its neighbour from mid point to mid point (see diagram). Measurements must be recorded in millimetres only.



Reading the tapes

Always read along the straight edge and in direction of arrows as indicated below.



Email or return forms to:

2- 4 Queen Street
Nunawading Victoria 3131, Australia
T (Int'l): +61 3 9915 8000
E: customer.service@jobskin.com
W: www.jobskin.com

Garment returns:

2- 4 Queen Street
Nunawading Victoria 3131, Australia

Go to www.jobskin.com
for measuring videos and tips

