

Patient details (Please print
 in black with CAPITAL letters)

Last name

First name

Date DD/MM/Y YYY



Calf Sleeve - 0015



Knee Length - 0101



Anklet - 0105

Interdigital Web
Spacer - 0041

Toe Piece - 1198B

Thong Toe Ankle/
Knee Length - 0539

Foot Glove - 0538

Measurements

Pleat			Pleat			
Left			Toe End	Right		
cm/mm				cm/mm		
			6			
			4 1/2			
			3			
			1 1/2			
			Heel 0			
			1 1/2			
			3			
			4 1/2			
			6			
			7 1/2			
			9			
			10 1/2			
			12			
			13 1/2			
			15			
Pleat			Knee End	Pleat		

Garment styles

Code	L	R	Description
0015			Calf Sleeve
0101			Knee Length
0105			Anklet (Garment length up to + 4 1/2 tape)
0041			Interdigital Web Spacer (to be worn over Foot Glove 0538)
1198B			Toe Piece ☐ Open Tips ☐ Closed Tips
0539			Thong Toe Anklet ☐ Open Tips ☐ Closed Tips
0539			Thong Toe Knee Length ☐ Open Tips ☐ Closed Tips
0538			Foot Glove ☐ Open Tips ☐ Closed Tips
0538			Individual Toe Knee Length ☐ Open Tips ☐ Closed Tips

Options and modifications

Code	L	R	Description
1169			Slant Inserts
0000			5cm Elastic
SBE			Silicon Backed Elastic ☐ 2.5cm ☐ 5cm
1190			Anterior Ankle Insert (powernet & lining)
1201			Seamless Sole ☐ Double Powernet ☐ Full Suede
1183			Suede Reinforcement ☐ Heel ☐ Ball
1159			Self Enclosed Toe
1160			Soft Enclosed Toe
1191			Slingback (for toe piece)
1164			Zippers (20cm or less)
1165			Zipper (20cm or more)
ZCT			Zipper Velcro® Closure Tab (top of zipper)
1163B			Zipper Fastening Aids (along length of zip)
1172A			Zipper Loop (leather)
0000			Zipper Lining ☐ Satin ☐ Elastic
0027			Padded Zipper Lining
1187			Lining-Heel
1189			Lining-Anterior Ankle
1188			Lining-Elsewhere (specify location) ☐ Satin ☐ Cotton
1196			Fully Lined Garment ☐ Soft Satin ☐ Cotton
			SilonTEX® lining – indicate location and size on either foot outline, TBSA map or printed photo
0027			Padding - Soft (indicate on TBSA Map) ☐ Pocket ☐ Attached
0027			Padding - Velfoam (indicate on TBSA Map) ☐ Pocket ☐ Attached
1202			Motif

Powernet colour selection guide

Fabric & Thread: Black, Brown, Beige, Light Beige, Purple, Navy, Royal, Light Blue, Red, Pink, White, Camouflage

Contrast Thread: Orange, Hot Pink, Lime, Yellow
Dual Stitching: nominate 2 colours

Fabric	Thread	Qty
1	1	
2	2	
3	3	
Must provide extra colour options		
1	1	
2	2	

NOTE: Powertex available in Beige and Black only

Comments

Leg zipper locations

	Location		Full Length (cm)	Or Specify Length
	L	R		
Lateral				
Medial				
Other				

Knee High, Anklet, Toe Piece

Measurement Form – Page 2 of 3

Patient details (Please print in black with CAPITAL letters)

Last name

First name

Date

DD/MM/YYYY



LEFT

☐ Open Tips
☐ Closed Tips

Desired toe lengths →

Length

Circumference of MTP (for toe piece only)

Circumference at end of toe piece

TOE CIRCUMFERENCES

Toe circumference 'B' should be 'base of toe' & 'A' should be 'at point of finished length' (mm)

(White boxes are lengths and grey boxes are circumferences)

RIGHT

☐ Open Tips
☐ Closed Tips

← Desired toe lengths

Length

Circumference of MTP (for toe piece only)

Circumference at end of toe piece

Note: FOOT OUTLINE (SEE: PG 3) MUST ACCOMPANY THIS FORM

Patient details (Please print in black with CAPITAL letters)Date

DD/MM/YYYY



Heel End

All measurements taken with no tension and in centimetres only.

Design options




Foot outline measuring tips

Centre the patient's foot on the foot outline page as marked. Ensure that the patient is standing or sitting with the foot flat on the ground. Mark the following clearly:


- Back of the heel
- Webspaces
- Toe lengths (if measuring a toepiece or glove)
- Where the garment will finish

Use the grey measuring tape for these garments. Measure in centimetres only.




Pleat instructions

When pleating the last tape, measure the distance between the pleated tape and its neighbour from mid point to mid point (see diagram). Measurements must be recorded in centimetres only.



Reading the tapes

Always read along the **straight edge** and in direction of arrows as indicated below.



Email or return forms to:

Unit 1, Level 1, 100 Station Street
Nunawading Victoria 3131, Australia
T (AU): 1800 772 577
T (NZ): 0800 765 627
T (Int'l): +61 3 9915 8000
E: customer.service@jobskin.com
W: www.jobskin.com

Garment returns:

Unit 1, Level 1, 100 Station Street
Nunawading Victoria 3131, Australia