

Thigh High, Waist High

Measurement Form - Page 1 of 3

Patient details (Please print in black with CAPITAL letters)



Jobskin

Working together for better solutions

+61 3 9915 8000

customer.service@jobskin.com

Jobskin.com

Last name	First name	Date
-----------	------------	------

Waist High, Two Legs, Closed Pubis
1101

Waist High, Capri / footless tights
1118

Waist High, Two Legs, Open Pubis
1102

Waist High, One Leg, Panty, Open Pubis
1112

Waist High, One Leg, Open Pubis
1103

Waist High, One Leg, Panty, Closed Pubis
1113

Bike Shorts, Closed Pubis
1119

Bike Shorts, Open Pubis
1122

Waist High Two Stumps
1133

Knee Sleeve
0015

Waist High, One Leg plus Stump, Closed Pubis
1134

Waist High, One Leg plus Stump, Open Pubis
1135

Waist High, One Stump,
1132

Chaps Style, Two Legs
0036

Chaps Style, One Leg
0035

Stump Support, (Above Knee)
1130

Panty Brief
1120

Stump Support (Below Knee)
1131

Thigh Length
0201

Thigh Sleeve
0019

Garment Styles				
Code	Qty	Description		
1101		Waist High, Two Legs, Closed Pubis		
1102		Waist High, Two Legs, Open Pubis		
1112		Waist High, One Leg, Panty, Open Pubis	☐ Left ☐ Right	
1103		Waist High, One Leg, Open Pubis	☐ Left ☐ Right	
1113		Waist High, One Leg, Panty, Closed Pubis	☐ Left ☐ Right	
1119		Bike shorts, Closed Pubis		
1122		Bike shorts, Open Pubis		
1133		Waist High, Two Stump		
0015		Knee Sleeve	☐ Left ☐ Right	
1134		Waist High, One Leg plus Stump, Closed Pubis	☐ Left ☐ Right	
1135		Waist High, One Leg plus Stump, Open Pubis	☐ Left ☐ Right	
1132		Waist High, One Stump	☐ Left ☐ Right	
0036		Chaps Style, Two Legs		
0035		Chaps Style, One Leg	☐ Left ☐ Right	
1130		Stump Support (Above Knee)	☐ Left ☐ Right	
1120		Panty Brief		
1131		Stump Support (Below Knee)	☐ Left ☐ Right	
0201		Thigh Length	☐ Left ☐ Right	
0019		Thigh Sleeve	☐ Left ☐ Right	
1118		Waist High, Capri Pants / Footless Tights		

Pubis Styles			
Male		Female	
1. Open	☐	1. Powernet & Lining	☐
2. Diagonal Fly	☐	2. Satin	☐
3. Gusset (fewer inseams)	☐	3. Mesh	☐
4. Pouch	☐	4. Open	☐
Options for male styles 2,3 & 4		Groin closure (for Body Brief only)	
Powernet	☐	Velcro® ☐ Snaplock ☐	
Powernet & Lining	☐		
Satin	☐		

Linings				
0040		Lining - Behind Knee	☐ Satin ☐ Cotton	
1177		Lining - Anterior Knee	☐ Satin ☐ Cotton	
1189		Lining - Anterior Ankle	☐ Satin ☐ Cotton	
1196		Fully Lined Garment	☐ Satin ☐ Cotton	
1188		Lining - Elsewhere (Indicate on TBSA Map)	☐ Satin ☐ Cotton	
		SilonTEX® lining - indicate location and size on either foot outline, TBSA map or printed photo		

Options and Modifications			
Code	L	R	Description
1177B			Knee Patch (Powernet)
1190			Anterior Ankle Insert (powernet & lining)
1201			Seamless Sole Standard for children up to 10 yrs old
1183			Suede Reinforcement ☐ Heel ☐ Ball
1159			Self Enclosed Toe
1160			Soft Enclosed Toe
1191			Stirrups
1169			Slant Inserts
1162			Suspenders - Attached ☐ 2.5cm ☐ 5cm
1107			Suspenders - Detached ☐ 2.5cm ☐ 5cm
SBE			Silicone Backed Elastic at Waist ☐ 2.5cm ☐ 5cm
0000			5cm Elastic Waistband
1163			Waistline Velcro® Attachment Tabs
1175			Double Abdominal Panel
1176			Uplift Panel
1161			Soft Abdominal Panel (Pregnancy)
1203			Postural Support (Indicate on TBSA Map)
0027			Padding - Soft (Indicate on TBSA Map) ☐ Pocket ☐ Attached
0027			Padding - Velfoam ☐ Pocket ☐ Attached
1202			Motif

Zippers			
1164			Waist Zippers
1165			Leg Zippers
ZCT			Zipper Velcro® Closure Tab
1163B			Zipper Fastening Aids (along length of zipper)
1172A			Zipper Loop (leather)
0000			Zipper Lining ☐ Satin ☐ Elastic
0027			Padded Zipper Lining

Comment

Thigh High, Waist High

Measurement Form - Page 2 of 3

Patient details (Please print in black with CAPITAL letters)



Working together for better solutions


+61 3 9915 8000

customer.service@jobskin.com

Jobskin.com

Last name	First name	Date
-----------	------------	------

Measuring points



Stump measurements

Medial to lateral	12				
Anterior to posterior	13				
Fold of Buttocks to Distal/Stump End, Right	11				
Fold of Buttocks to Distal/Stump End, Left	11				

Body Measurements

	Circumference (mm)				Heights (mm)			
Desired top of Support	1				7			
Waist	2				8			
Midpoint between 2 and 4	3							
Largest Part of Buttocks	4							
Gluteal Fold, Left	5							
Gluteal Fold, Right	6							
Fold of Buttocks to Distal End					9			
Fold of Buttocks to Top of Support					10			

Personalisation

Fabric Selection

- ☐ Conformtex
☐ Powernet
☐ Powertex

Colour Selection

Refer to Jobskin swatch or [website](#) for fabric and thread colour choices

Fabric Colour	Thread Colour	Qty
1	1	
2	2	
3	3	

Lower limb circumferences

Indicate tape which is closest to popliteal crease

Pleat	Pleat
Left	Right
To end	To end
mm	mm

			7 1/2		
			6		
			4 1/2		
			3		
			1 1/2		
			Heel 0		
			1 1/2		
			3		
			4 1/2		
			6		
			7 1/2		
			9		
			10 1/2		
			12		
			13 1/2		
			15		
			16 1/2		
			18		
			19 1/2		
			21		
			22 1/2		
			24		
			25 1/2		
			27		
			28 1/2		
			30		
			31 1/2		
			33		
			34 1/2		

Pleat

Pleat

Zipper locations

	Location		Full Length	Or Specify Length (cm)
	L	R		
Lateral Aspect				
Medial Aspect				
Other:				
Waist Zippers				

Comments

Thigh High, Waist High

Measurement Form - Page 3 of 3

Patient details (Please print in black with CAPITAL letters)

Last name

First name

Date



Jobskin
Working together for better solutions
+61 3 9915 8000
customer.service@jobskin.com
Jobskin.com

2726252423222120191817161514131211109876543210

cm

0

1

2

3

4

5

6

7

8

9

10

11

12

13

14

15

16

16

For open tips - indicate on foot outline desired length

Please ensure all bases of webspaces are clearly marked

Place third toe on line

Toe end

Heel end

*Velcro is a registered trademark of VELCRO BVBA

All measurements taken with no tension and in millimetres only

Design options



Bike Short



Panty Brief



Waist-high
2 x leg



Waist-high
1 x leg



Waist-high
1 x leg w panty




Thigh-high



Chap-style



Gusset



Open Crotch



Childs Open Crotch



Diagonal Fly



Pouch



Ostomy Opening
(not recommended)



Uplift Panel



Postural Support
(many variations)




Pregnancy Panel



Suspenders

Leg measuring tips


Use blue JOBSKIN tape for full length lower limb measurements



Use the pink JOBSKIN tape for paediatrics


Pleat instructions

When pleating the last tape, measure the distance between the pleated tape and its neighbour from mid point to mid point (see diagram).



Reading the tapes

Always read along the straight edge and in direction of arrows as indicated below.



Email or return forms to:

2- 4 Queen Street
Nunawading Victoria 3131, Australia
T (Int'l): +61 3 9915 8000
E: customer.service@jobskin.com
W: www.jobskin.com

Garment returns:

2- 4 Queen Street
Nunawading Victoria 3131, Australia

Go to www.jobskin.com
for measuring videos and tips

