

Patient details (Please print in black with CAPITAL letters)

Working together for better solutions

+61 3 9915 8000
customer.service@jobskin.com
Jobskin.com

Last name


First name

Date


Measuring points

Measurements


Circumferences



Lengths: Back




Lengths: Front



Vertical Linear Lengths Horizontal Contour Lengths

Protruding abdomen



Vertical Linear Lengths Horizontal Contour Lengths

Protruding abdomen

Neck circumference	1				
Left shoulder circumference	2				
Right shoulder circumference	3				
Chest circumference	4				
Circumference over nipple line	5				
Circumference just under breast (level of breast attachment to chest)	6				
Waist circumference	7				
Circumference over largest part of belly	8				
Circumference at end of support	9				
C7 to nipple line	10				
C7 to just under breast	11				
C7 to waist	12				
C7 to largest part of belly	13				
C7 to end of support	14				
Side seam to side seam under axilla on back	15				
Side seam to side seam at waist on back	16				
Shoulder width	17				
Singlet/Sleeveless, garment shoulder width	17a				
Top arm tape to neck (arm down by side)	18				
Under breast to waist on front	19				

Contour, sternal notch to end of support	20				
--	----	--	--	--	--


















Patient to wear well fitted bra during measurement of Breast Support.

Patient details (Please print in black with CAPITAL letters)

Last name	First name	Date
-----------	------------	------

Upper limb circumferences



mm				mm	
					
					
					
					
					
					
					
					
					
					
					
					
					
					
					
					
					



Note: Only pleat tapes at ends ie:at wrist and/ or axilla

Arm zipper locations

	Location		Full Length	Or Specify Length (cm)
	L	R		
Lateral/Radial Aspect				
Medial/Ulnar Aspect				
Other:				


Personalisation

Fabric Selection <input type="checkbox"/> Confortex <input type="checkbox"/> Powernet		Colour Selection Refer to Jobskin swatch or website for fabric and thread colour choices	
Fabric Colour	Thread Colour	Qty	
1	1		
2	2		
3	3		

Comments

***Patient photos assist with garment design.
 For complex, unusual or very asymmetrical breasts, please use separate Cutom Bra Cup Form.***

Design options



Design choice definition

Funnel neck: an extension of the vest up the neck without the circumferential neckline seam. Measure height as per attached powernet turtleneck.

Full neck scoop: This option will increase the width and depth of the neckline as well as decreasing compression over these areas.

Shoulder bracing: A diagonal pull to help extend the shoulders backwards and downwards to assist in improving posture as well as garment fit over the anterior chest wall.

Velcro® Closure: This closure option allows adjustability and is often used for children. It cannot be used as a front closure for women when bra cups are included.

Torso measuring tips

Neck circumference (1): Measured midway between the head and the shoulders.

Round neck back: Measure from C7 to depth of scoop.

Round neck front: Measure from sternal notch to depth of scoop.

NB. Round neck. This option is not recommended when treating the upper torso as a large scoop will significantly reduce the pressure around the shoulder, upper chest and back.

Attached powernet turtleneck: Measure height from sternal notch.

Shoulder circumference (2 & 3): With patient standing with hand on hip, place tape measure on top of shoulder directly above axilla, and pass through under the arm. Then lower arm down and continue tape to shoulder. You will notice the measure 'ease'.

Arm circumference: Using Purple JOBSKIN arm tape. or Pink JOBSKIN tape for paediatrics.

For **Elbow contracture seam**, record the degrees of **deficit** in range of motion

Chest circumference (4): Measured at the level of the axilla with arms placed down to the side.

Nipple line (5): For women please provide bra cup size. Circumference over nipple line should be measured with patient wearing a bra.

Side seam to side seam (15 & 16): A horizontal contoured measure originating directly below the midpoint of one axilla across the back to the midpoint under the other axilla (15) is positioned at the height of chest circumference (4), whilst (16) is positioned at the waist.

Nappy straps: Measure over the nappy from waist at back to waist at front.

Vest length (14 & 20): We recommend the vest extend well beyond the waist for anchorage to help prevent garment riding up.

Trachea/Ostomy opening: If opening is required, it is recommended to fit the garment first, mark the position required, then return to the factory for finishing. NOTE: Ostomy openings are not recommended because of the risk of herniation. Obtain medical clearance.

Email or return forms to:

2- 4 Queen Street
Nunawading Victoria 3131, Australia
T (Int'l): +61 3 9915 8000
E: customer.service@jobskin.com
W: www.jobskin.com

Garment returns:

2- 4 Queen Street
Nunawading Victoria 3131, Australia

Go to www.jobskin.com
for measuring videos and tips

