

Patient details (Please print in black with CAPITAL letters)




Working together for better solutions

+61 3 9915 8000
customer.service@jobskin.com
Jobskin.com

Last name	First name	Date
-----------	------------	------

How to complete hand outline page

CLEARLY MARK END OF GARMENT ON FINGERS



Gauntlet - MCPs to Wrist **0505**



Glove to Wrist **0535**



Interdigital Web Spacer **0536**



Mitten (for amputation at MCP level) **0537**



Gauntlet to Axilla **0502**



Glove to Axilla **0532**



Gauntlet to Elbow **0516**



Glove to Elbow **0534**



Gauntlet to Axilla with Shoulder Attachment **0504**

Garment styles

Code	L	R	Description
0505			Gauntlet - MCPs to wrist
0535			Glove to Wrist <input type="checkbox"/> Open Tips <input type="checkbox"/> Closed Tips
0536			Interdigital Web Spacer <input type="checkbox"/> Thumb <input type="checkbox"/> No Thumb
0537			Mitten to Wrist (amputee)
0502			Gauntlet to Axilla
0532			Glove to Axilla <input type="checkbox"/> Open Tips <input type="checkbox"/> Closed Tips
0516			Gauntlet to Elbow
0534			Glove to Elbow <input type="checkbox"/> Open Tips <input type="checkbox"/> Closed Tips
0504			Gauntlet with Shoulder Attachment
0025			Individual Finger Stall <input type="checkbox"/> Open Tips <input type="checkbox"/> Closed Tips <input type="checkbox"/> Wrist Band

Option and modifications

Code	L	R	Description
			Seam-free Thenar Eminence
0021			Reinforced Palm (Powernet)
1183			Palm Grip <input type="checkbox"/> Silicone <input type="checkbox"/> Suede
1188			Positional Reinforcement (indicate on hand outline)
1169			Slant Inserts
0027			Padding - Soft <input type="checkbox"/> Pocket <input type="checkbox"/> Attached <small>(indicate on TBSA Map or hand trace)</small>
0027			Padding - Velfoam <input type="checkbox"/> Pocket <input type="checkbox"/> Attached <small>(indicate on TBSA Map or hand trace)</small>
1202			Motif
SBE			Silicone Backed Elastic
1190			Anterior Elbow Insert
0000			Left Elbow Contracture Seamdegrees
0000			Right Elbow Contracture Seamdegrees
GCCINSERT			Contrast Colour Inserts (Specify 3 choices) 1..... 2..... 3.....
0021			Lining-Palm
1167			Lining-Inside Elbow <input type="checkbox"/> Satin <input type="checkbox"/> Cotton
1168			Lining-Full Elbow <input type="checkbox"/> Satin <input type="checkbox"/> Cotton
1167B			Lining-Posterior Elbow <input type="checkbox"/> Satin <input type="checkbox"/> Cotton
1188			Lining (specify location) <input type="checkbox"/> Satin <input type="checkbox"/> Cotton
1197			Fully Lined Garment <input type="checkbox"/> Soft Satin <input type="checkbox"/> Cotton
			SilonTEX® lining - Indicate location and size on either hand outline, TBSA map or printed photo
1164			Zipper
1172A			Zipper Loop (leather)
ZCT			Zipper Closure Velcro® Tab (top of zipper)
1163B			Zipper Fastening Aids (along length of zipper)
0000			Zipper Lining <input type="checkbox"/> Satin <input type="checkbox"/> Elastic
0027			Padded Zipper Lining
1188			Smart Touch Finger Pad (specify location in comments)

Personalisation

Fabric Selection	Colour Selection
<input type="checkbox"/> Conformtex	Refer to Jobskin swatch or website for fabric and thread colour choices
<input type="checkbox"/> Powernet	

Fabric Colour	Thread Colour	Qty
1	1	
2	2	
3	3	

Hand/arm zipper locations

	Location		Full Length	Or Specify Length (cm)
	L	R		
Dorsal				
Ulnar				
Palmar				
Other :				

Comments

Patient photos assist in garment design.

Glove, Gauntlet


Measurement Form - Page 2 of 3

Patient details (Please print in black with CAPITAL letters)

Last name

First name

Date




Arm circumferences

Pleat	Wrist end	Pleat
Left (mm)		Right (mm)
	-3	
	-1 1/2	
	0	
	1 1/2	
	3	
	4 1/2	
	6	
	7 1/2	
	Elbow 9	
	10 1/2	
	12	
	13 1/2	
	15	
	16 1/2	
	18	
	19 1/2	
Pleat	Axilla End	Pleat

Note: Only pleat tapes at end ie: at wrist and/or axills

Shoulder attachment measurements			
1. Distance last arm tape to inside edge of bra strap	2. Top of shoulder, circumferentially under opposite breast	3. Top of shoulder, circumferentially under opposite axilla	4. Top of shoulder, diagonal to waist
Bra Loop (mm)	Female Body Strap (mm)	Male/Child's Body Strap (mm)	Waist Attachment (mm)
1	2	3	4
	2a	3a	4a




Glove, Gauntlet

Measurement Form - Page 3 of 3

PLEASE ENSURE THIS PAGE PRINTS TO CORRECT SCALE BEFORE TAKING HAND TRACE

Patient details (Please print in black with CAPITAL letters)



All measurements taken with no tension and in millimetres only

Design options



Glove



Gauntlet



Interdigital
Webspacer



Finger Stall
wrist band



Anterior
Elbow insert



Reinforced
Palm



Palm Grip
(silicone)



Palm Grip
(suede)



Slant inserts



Thenar Eminence
Free Thumb

Gloves & gauntlets definition, hints and tips

Gauntlet measurements required: palm, thumb, wrist and any arm measurements needed. The hand outline must be completed.

The important points to mark on the hand trace are the web spaces, wrist and finishing point of glove.

If fingers are contracted, measure from the web space along the side of the fingers and follow the curve to the desired length of each finger.


A glove to wrist may extend up to tape 3.

Finger Circumferences: Take with a JOBSKIN finger tape and measure in skin contact - not tightly.

Finger Lengths: Can be taken with the JOBSKIN Handy-card.

Open Tips: Measure only the length desired from the Finger Web Space (palmar aspect).

Slant Inserts: Enable a snug fit on the dorsal surface and give extra pressure in the web spaces.



Feed measurement through slot at top of tape


Use the finger tape to measure finger circumferences of hand. Use the purple tape for the arm.

Web Spacers: Are designed to give extra pressure in the web spaces and are worn over gloves. They can be ordered with or without a thumb piece.

Zippers: May be used to allow for inserts, e.g. silicone gel sheets.

Arm measuring tips

With the arm straight and horizontal, proceed to measuring using a purple tape.



Use pink JOBSKIN tape for paediatrics.


For a Full Length Sleeve begin by securing tape 9 directly around the elbow. The last tape should fall just below the wrist bone (distal to styloid process). Pleating at the spine occurs ONLY at the wrist and axilla if required.

For a Forearm Sleeve begin tape at wrist and pleat at proximal end if required.

For Elbow Contracture Seam record the degrees of deficit in range of motion.


Pleat instructions

When pleating the last tape, measure the distance between the pleated tape and its neighbour from mid point to mid point (see diagram).



Reading the tapes

Always read along the straight edge and in direction of arrows as indicated below.



Email or return forms to:

2- 4 Queen Street
Nunawading Victoria 3131, Australia
T (Int'l): +61 3 9915 8000
E: customer.service@jobskin.com
W: www.jobskin.com

Garment returns:

2- 4 Queen Street
Nunawading Victoria 3131, Australia

Go to www.jobskin.com
for measuring videos and tips

